



Ultrasonic Nebulizer

PRIZMA _ Germany
Technical Data

Model : PROFI SONIC

- Nebulization rate: 0.0-3.0 ml/min
- Air volume: max 20 l/min
- Particle size: 0.5-5.0 microns
- Capacity of medication cup: 5-150 ml
- Capacity of water tank: 300 ml
- Low drug residue: 1 ml
- Ultrasonic frequency: 1.7 MHz
- Noise level: <35 dB (A)
- Operating voltage/frequency: 100-240 V AC, 50/60 Hz
- Power consumption: 55 W
- Dimensions: 260 x 200 x 250 mm
- Weight: 1.8 kg

CONTENTS OF THE KIT:

PROFI SONIC nebulizer, inhalation hose (75 cm), mouthpiece, 3 medication cups (one already installed), adult mask, child mask, bacterial filter, power cord and instruction manual.

Optional accessories available for different types of inhalation therapy.

Touch screen & phone
for parameters setting

Large number of patients
continuous inhalations

In every position
nebulization



Specifications

Technical data

This device fulfils the provisions of the EC directive 93/42/EEC (Medical Device Directive).

Product name	PRIZMA Ultrasonic Nebuliser
Model	PROFI SONIC
Product class	Ila
Power source	100-240V~, 50-60Hz
Power consumption	55W
Fuse	2 x T2AH/250VAC
Electrical safety	I, Type B
Electromagnetic emissions	In accordance with EN 60601-1-2 standard
Noise level	<35dBA
Ultrasonic frequency	1.7MHz
Nebulisation rate*	0-3 ml/min, adjustable
Particle size	0.5-5 µm
Air volume	Maximum 20 l/min, adjustable
Amount of cooling water	300 ml
Capacity of medication cup	150 ml (min. 5 ml)
Clock battery	CR2032, 3V
External dimensions	260 x 250 x 200
Weight of the main unit	1.8kg
Operating conditions	15 - 40°C (59 - 104°F), 30 - 85% RH, 700-1060 hPa
Transport/storage conditions	0 - 45°C (32 - 113°F), 30 - 85% RH, 700-1060 hPa
Package contents	Main unit, silicone inhalation hose (90 cm), mouthpiece set, 3 medication cups (one already installed), 1 child mask, 1 adult mask, 1 bacterial filter (already installed), power cord, Instruction Manual (with warranty card)

* May change depending on the type of medication

- Specifications and appearance is subject to change without prior notice.

Symbols:



= Electrical safety: Type B



= Read the instruction manual carefully



= Correct disposal of this product (Waste electrical & electronic equipment)

Manufacturer:
PRIZMA D.O.O.
Kumanovska 8
34000 Kragujevac, Serbia



CE 1304

tel: 07711011137-07715553335
Email: info@almanara-iq.com
website : www.almanara-iq.com

بغداد / باب المعظم / خلف جامع عادلة خاتون

المناصرة
مكتب المناصرة الزرقاء العلمي