Genex Laboratories

Limited liability Company 300 Avenue south, Florida, USA



THROMBO

Semi Automated Coagulation analyzer

Technical Specifications

Wavelength: LED 470nm

Memory: 10,000 results storage

Input: Keypad

Sample volume: 20-40µL

Reagent volume: 20-40µL

Electronic pipette: 20-200µL

Available dimensions: Seconds, mg/dL, g/L, U/L

Incubation temperature: 37.0±1.0⊠

Relative humidity: ≤85%

Precision: PT, APTT, FIB \leq 5%; TT \leq 8% (Normal sample)

Deviation between channels: ≤ 5%

Power supply: AC 100-240V 50/60Hz

2 independent channels

Display: 5-inch screen

Cuvette incubation positions: 12

Reagent incubation positions: 5

Size: 345mm×325mm×145mm

Net weight: 4.4kg

Tel: 07736768000

Email: sales@almanara-iq.com website: www.almanara-iq.com

بغداد / باب المعظم / خلف جامع عادلة خاتون



